

## **ApneaLink**<sup>™</sup>Air

Home Sleep Testing Device



## Home sleep testing made easy

ApneaLink<sup>™</sup> Air provides the performance and reliability of our trusted ApneaLink products in a sleeker, easier to use device. With its one-touch operation and intuitive Start/Stop button, it's never been easier for patients to set up and use successfully. And the test complete indicator helps minimize wasteful retesting and increase patient satisfaction.



# ApneaLink Air: Our most advanced home sleep testing solution

With its stylish design, ApneaLink Air provides performance and reliability in a compact, lightweight and easy-to-use home sleep testing (HST) device. A low-cost, type III HST device, the ApneaLink Air is capable of recording up to five channels of information: respiratory effort, pulse, oxygen saturation, nasal flow and snoring. It also features longer recording time and increased storage capacity over previous ApneaLink products.



### Access to detailed signals overnight

ApneaLink Air software provides clinicians access to a more in-depth view of their patients' recordings.



## ApneaLink Air benefits

- **Improved usability** for both patients and healthcare professionals with:
  - One-touch operation with an intuitive Start/Stop button
  - Longer recording time (~12 hours) and increased storage capacity (~48 hours)
  - Test complete light provides feedback on recording time to users
  - Good signal indicators to help successful set ups

#### • Reliable automatic analysis

- Includes apnea–hypopnea index (AHI), hypopnea index (HI), flow limitation, snoring and oxygen desaturation index (ODI)
- Manual scoring possibility for more detailed patient data
- Configurable analysis parameters allow adjustment of thresholds

- Validated results meet AASM and CMS definitions for hypopnea scoring guidelines
- Cheyne–Stokes probability detection helps determine when to refer patients for further in-lab diagnosis
- Differentiation between obstructive and central apneas
- **Printable patient instructions** when programming the device
- Flexible report options
  - **Email summary and signal reports** to send to referral physicians or other relevant parties
  - Extended report contains additional overview of respiratory data
  - **Prescription page** streamlines process for healthcare professionals

## ApneaLink Air Complete Kit



## Comprehensive reporting

ApneaLink Air automatically analyzes and derives AHI, flow limitation and snoring and later automatically generates a simple, easy-to-interpret report with a color-keyed AHI or Risk Indicator for the clinician to review.

First name:       Cheyne Stokes       Patient ID: DOB:       D/1/2/1958         Street:       170.18 cm       000         Street:       170.18 cm       000         Street:       170.19 cm       000         Street:       170.19 cm       000         Street:       170.19 cm       000         Street:       21:00       170.19 cm         Phone:       21:13       170.19 cm         Date:       21:13       170.19 cm         Start:       21:23       170.19 cm         Date:       21:00       170.19 cm         Duration:       8 h 43 min       100 cm         Duration:       9 h 43 min       100 cm         AHI*       100 cm       100 cm         Result       100 cm       100 cm         Athi:       22.3       5/n         Ath:       23.3       5       Breaths:       6857         Apnea index:       20.4       5/n       Average breaths per minute [bpm]:       14.45         All:       0       0       0       00 cottructure apneas:       116 (72%)         All:       14.9       Central apneas:       116 (72%)       117 (15)         Aporead, 6 5/n	Treating physician			Referral to	
Date:       21/01/2009         Start:       21/13.         End:       21:3.         Duration:       8 h 43 min         Start:       20:564.         Duration:       8 h 43 min         Object       Duration:         AHI*       Image: Ima				DOB: Height: Weight:	170.18 cm 92.41 kg
Normal range       Suspected pathological breathing disorder         Result (33)         * see Circical Quide for abbreviations and Realified standard parameters         Analysis (Flow evaluation period: 7 h 55 min / SpO2 evaluation period: 8 h 32 min) indices       Normail       Result         Analysis (Flow evaluation period: 7 h 55 min / SpO2 evaluation period: 8 h 32 min) indices       Normail       Result         AH1::       32.9       < 5 / h       Average breaths per minute [bpm]:       14.45         AH1::       0.6       < 5 / h       Apreasi:       0 (0%)         UAI:       0       Unclassified apneas:       0 (0%)         OAI:       5.6       Obstructive apneasa:       118 (72%)         MAI:       0.1       Mixed apneas:       118 (72%)         MAI:       0.1       Flow lim. Br. with sn (FL):       797         % Flow lim. Br. with Sn (FS):       0       < Approx.40       Flow lim. Rr. with Sn (FL):       33         Shoring events:       66       ODI Oxygen Desaturation:       236       5 / h       No. of desaturations:       2	Date: 21/01/20 Start: 21:13 . End: 05:56 .			Start: End:	05:54 .
Result (33)         'see Clinical Guide for abbreviations and ResMed standard parameters         Analysis (Flow evaluation period: 7 h 55 min / SpO2 evaluation period: 8 h 32 min)         Indices         Normal Result         All ': 32.9       5 / h       Average breaths per minute [bpm]:       14.45         All ':       33.8       < 5       Breaths:       6857         Apnea index:       20.6       < 5 / h       Apneas:       0 (0%)         OAI:       5.6       Obstructive apneas:       4 (27%)         CAI:       14.9       Central apneas:       118 (72%)         MAI:       0.1       Mixed apneas:       1 (1%)         Hypopnea index:       12.3       < 5 / h       Hypopneas:       97         % Flow lim. Br. with out Sn (FL):       12       < Approx.40       Flow lim. Br. without Sn (FL):       792         % Flow lim. Br. with Sn (FS):       0       < Approx.40       Flow lim. Br. with Sn (FS):       3         Cobit Oxygen Desaturation:       72       -       Saturation <= 90% :       234 min (45%)         Lowest desaturation:       72       -       Saturation <= 85% :       77 min (15%)         Dol Oxygen Desaturation:       72       <			A		
Analysis (Flow evaluation period: 7 h 55 min / SpOz evaluation period: 8 h 32 min) Indices       Normal       Result         AHI":       32.9       < 5 / h       Average breaths per minute [bpm]:       14.45         AHI":       33.8       < 5       Breaths:       6857         Apnea index:       20.6       < 5 / h       Apneas:       163         UAI:       0       Unclassified apneas:       0 (0%)         OAI:       5.6       Obstructive apneas:       44 (27%)         CAI:       14.9       Central apneas:       118 (72%)         MAI:       0.1       Mixed apneas:       1 (1%)         Hypopnea index:       12.3       < 5 / h       Hypopneas:       97         % Flow lim. Br. without Sn (FL):       12       < Approx. 60       Flow lim. Br. without Sn (FL):       792         % Flow lim. Br. with Sn (FS):       0       < Approx. 40       Flow lim. Br. with Sn (FS):       3         ODI Oxygen Desaturation Index*:       33.5       < 5 / h       No. of desaturations:       286         Average saturation:       91       94% - 98%       Saturation <= 85% :       77 min (15%)         Lowest desaturation:       72       90% - 98%       Saturation <= 80% :       7 min (1%)         Baseline Satura		••••		•••••	
Analysis status: Analyzed automatically	Indices AHI*: RI*: Apnea index: UAI:	32.9 33.8 20.6 0 5.6	Normal < 5 / h < 5	Result Average breaths per minute [bp Breaths: Apneas: Unclassified apneas: Obstructive apneas: Central apneas:	6857 163 0 (0%) 44 (27%) 118 (72%)
Aprea [20%; 10s; 80s; 1.0s; 20%; 80%; 8%]; Hypopnea [70%; 10s; 100s; 1.0s]; Snoring [8.0%; 0.3s; 3.5s; 0.5s]; Desaturation [4.0%]; CSR [0.50]	CAI: MAI: Hypopnea index: % Flow lim. Br. without Sn (FL): % Flow lim. Br. with Sn (FS): ODI Oxygen Desaturation Index*: Average saturation: Lowest desaturation: Lowest saturation: Baseline Saturation: Minimum pulse: Maximum pulse: Average pulse:	12.3 12 0 33.5 91 72 98 59 71 62	< Approx. 60 < Approx. 40 < 5 / h 94% - 98% - 90% - 98% % > 40 bpm < 90 bpm bpm	Flow lim. Br. without Sn (FL): Flow lim. Br. with Sn (FS): Snoring events: No. of desaturations: Saturation <= 90% : Saturation <= 85% : Saturation <= 80% :	97 792 3 66 286 234 min (46%) 77 min (15%)



Home Sleep Testing Device

# Ordering information and product codes

### Accessories

### ApneaLink Air complete kit

- Software installation CD
- ApneaLink Air patient instructions
- USB download cable
- 1 reusable belt
- Oximeter reusable soft sensor (8000SM) Medium
- XPOD LP oximeter
- XPOD LP oximeter clip
- 3 nasal cannulas
- Carry bag
- 2 AAA batteries

### **Oximeter components**

XPOD LP oximeter	22374
XPOD LP oximeter clip	22371
NONIN 6000CA single use oximetry sensors – pk 24	22337
NONIN oximeter reusable soft sensor (8000SS) – Small	70567
Small	/050/
NONIN oximeter reusable soft sensor (8000SM) – Medium	70568
NONIN oximeter reusable soft sensor (8000SL) –	
Large	70413

22354

### **Nasal Cannulas**

Nasal cannula – pk 25	70388
-----------------------	-------

EasySense <sup>™</sup> Respiratory	
Effort Sensor	

22321

#### Other accessories

ApneaLink Air reusable belt	629052
ApneaLink Air carry bag	22373
ApneaLink Air software CD V10.1	22372
USB cable for ApneaLink Air	22375



For more information: Accessories Guide

Please contact your ResMed representative for more information.



## **ApneaLink**<sup>™</sup>Air

Home Sleep Testing Device

## ApneaLink Air Technical Specifications

Recorder	Download application
Enhanced hardware	PC download
EasySense respiratory effort sensor	Internal battery
Oximetry	AAA alkaline batteries
Light indicators	Internal memory
- Test complete indicator	Recording period: 48 hours
<ul> <li>Respiratory flow indicator</li> <li>Effort sensor connection indicator</li> <li>Oximetry connection indicator</li> </ul>	Internal clock
	Dimensions
Signal Recording	Recorder: 62 x 102 x 30 mm (2.4" x 4" x 1.2")
Respiratory effort	Pulse oximeter: 53 x 20 x 15 mm (2.1" x 0.8" x 0.6")
Respiratory flow	Recorder weight: 66 g (2.3 oz)
Snore	
Blood oxygen saturation	
Pulse	

Battery voltage

