


QUICK SETUP GUIDE

SECULIFE BP_{PRO}



SECULIFE BP_{PRO} is a microprocessor based, high precision, Non-Invasive Blood Pressure (NIBP) simulator. The unit is small, easy to use and has multiple features to fit many different applications. The SECULIFE BP_{PRO} offers invasive blood pressure, temperature, arrhythmias and a leak rate test mode. The graphical display provides multiple screens containing pressure in mmHg, a plot of the overall pressure, or a close-up of the BP waveforms. SECULIFE BP_{PRO} features 9 basic test modes.

SECULIFE BP_{PRO}



Fingersim SET



The FingerSim pulsoximeter-testing system enables the user to test pulsoximeters and sensors with three simulated light absorption conditions. These absorption conditions imitate typical oxygen saturations.

Three FingerSims:

SpO₂ 80 %

SpO₂ 90 %

SpO₂ 97 %

Incl. case, operating instructions and SECULIFE spare holder.

SECULIFE Pulse Oximeter Module



Finger holder

Simulates heart frequency in combination with the basic device

Running a test

The SECULIFE BP_{PRO} includes 9 Basic Test modes and 3 Sub-Test Modes. This section will guide you through each of the tests and their basic operations.

The Main Tests are accessible with a single key. The  key will scroll through the following tests in a continuous loop:

- Adult 120/80
- Adult 120/80 w/ Pace
- Adult High 190/120
- Adult Low 80/40
- Neonatal 70/40
 - Alarm Test
- Arrhythmia Sequence
 - Leak Test
 - Manometer

Basis Test Modes

NIBP:

The first five test modes deal with various NIBP setups. To run an NIBP simulation, the cuff and monitor are connected to the pressure input. Then the measurement is initiated by the monitor and the SECULIFE BP_{PRO} unit will show the proper waveform based on the cuff pressure provided by the monitor and the selected simulation.


The NIBP output mode can be changed by pressing the  key. Once the desired operating mode is selected, the output will begin when the correct pressure is detected automatically.

Test-Modes


There are 5 selectable NIBP Basic Test Modes: adult, adult w/ pace, adult high, adult low and neonatal. The following displays show examples for the selectable modes:

NOTE: The screens for the first test will have the individual components labeled. The component labels for subsequent tests are the same.


Adult 120/80 with Pace



Adult High 120/80



Adult Low 80/40



Neonatal 70/40



Examination of blood pressure simulator

- 1.) Connection of SECULIFE BP_{PRO} and DUT (monitor) the SECULIFE BP_{PRO} blood simulator is connected between the blood pressure sleeve and the DUT monitor. A T-piece is inserted in order to scan the pressure.



Detail: magnified T-piece

- 2.) Functional check – blood pressure

In order to test the correct function, the blood pressure values set at the SECULIFE BP_{PRO} (here: 120/80) are compared with the display values of the DUT (here: 123/80).



Detail: SECULIFE BP_{PRO} display



Detail: DUT display (monitor)

Examination of invasive blood pressure (IBP)

The output circuit is fully isolated and capable of switching between the two standard sensitivities (5 μ V/V/mmHg and 40 μ V/V/mmHg).

1.) Target values

The same specifications apply as for NIBP simulation.

2.) Settings

At the SECULIFE BP_{PRO} Display, it is possible to set the transformer ratio as well as the offset for the monitor.

3.) Before starting the test, some patient monitors may require a zero-balance with the SECULIFE BP_{PRO}.

The IBP simulation will start automatically, as soon as both instruments will be switched back to simulation- respectively to measuring mode.



Detail: SECULIFE BP_{PRO} display

Temperature test


- 1.) Connection of temperature test cable and SECULIFE BP_{PRO} for measuring the temperature, a standard temperature test cable is connected between SECULIFE BP_{PRO} and the DUT.

The temperature setting can be selected in the **SETUP** mode.




- 2.) Function test

The temperature at the SECULIFE BP_{PRO} (here: 35 °C) is compared with that of the DUT (here: 35.1 °C) in order to check the function of the temperature display.



Detail:: DUT (monitor) display



Detail: SECULIFE BP_{PRO} display

Examination of SpO₂ (oxygen saturation)

1.) Connection for SpO₂ measurement

A connection is established between SECULIFE BP_{PRO} and the holder via a 7 pin mini-DIN cable.



2.) Checking the SpO₂ function with the "SECULIFE pulse oximeter module". The nominal value of the FingerSim™ (here: 97 %) should approximately correspond to the value displayed at the DUT (here: 98 %).

NOTE: Please observe the calibration date of the FingerSim™.



Detail: DUT display (monitor)







Figure: magnified FingerSim 97 % SpO₂

SpO₂ temperature dependence

Changes in the ambient temperature have an impact on the light absorption characteristics of the FingerSim which results in slight temperature induced changes in the simulated SpO₂ values. Each FingerSim has been calibrated at 22.5 °C. If the ambient temperature ranges between 19.7 °C and 25.3 °C, it is not necessary to modify the simulation to be expected. If the ambient temperature is above or below this range, however, the diagrams shown on this page should be used to modify the simulated SpO₂ value to be expected.




Leakage test

1.) Connection

A connection is established between SECULIFE BP_{PRO} and the pressure sleeve by means of the T-piece. Moreover, for conducting the leakage test, it is necessary to fit the T-piece with a pressure ball. The pressure ball, including discharge valve, is connected to the T-piece by means of a Luer Lock tube to this end.

Note:

The pressure may not exceed **500 mmHg!**



2.) Settings

With the **SELECT KEY** you can scroll through the test modes and select the leakage test mode. This function provides a standardized pressure leakage test. The pressure input is connected with the system to be monitored. After that, the system is put under pressure. The leakage test is started by pressing the **RESET KEY**. The device locates the pressure drop and indicates it. Additionally, the time spent since the beginning of the test is indicated. Furthermore, the leakage rate in mmHg per minute is calculated and displayed.



Examination of ECG simulation

1.) Connection of SECULIFE BP_{PRO} with the DUT


The SECULIFE BP_{PRO} and the DUT (monitor) are connected with each other by means of a patient cable including snap lock clips. Take care of the color coding on the housing.



2.) Settings and Examination

Use the **SELECT** key to cycle through the test modes until "arrhythmia sequence" appears. The purpose of this function is to provide a timed sequence of some of the more common arrhythmias.

The 6 selectable heart rhythm disturbances (arrhythmias) of the SECULIFE BP_{PRO} are shown on the DUT display for verification purposes.



Examination of respiration


1.) Connection

Establish a connection between SECULIFE BP_{PRO} and the DUT (monitoring) for examining the respiration. The respiration is only issued via the LA-ECG cable. The connection is identical with the ECG simulation.



2.) Function test


To test the function, the value (here: 35 brpm) shown on the display of the DUT (here: 35 brpm) is compared to the value issued by the SECULIFE BP_{PRO} (here: 35 brpm).



Excerpt: magnified display

General view of overall setup

The entire function tests described before can be performed simultaneously.



GMC INSTRUMENTS



GMC-I Messtechnik GmbH

Südwestpark 15 □ 90449 Nürnberg □ Germany

TEL +49 911 8602-111 □ FAX +49 911 8602-777

www.gossenmetrawatt.com □ info@gossenmetrawatt.com