

### THE MOST EFFECTIVE NON-INVASIVE APPROACH TO COOLING THERAPY

#### The Optimal Solution for Therapeutic Hypothermia

Cure Wrap is an integral part of the CritiCool system. The Wrap covers the Patient in a suit-like matter, fastened with Hook & Loop strips. The Wrap can be opened easily for Patient Care and then refastened without impairing its use. Water Channels in the Wrap grant a constant flow of water to facilitate energy transfer to the Patient. The Water Temperature is regulated according to feedback from Core and surface Temperature Sensors placed on the Patient.

#### CureWrap - Step by Step Application



# CUREWRAP TYPES

**CUREWRAP GUARANTEES CLOSE, SAFE CONTACT WITH YOUR PATIENT'S SKIN DURING THERAPEUTIC HYPOTHERMIA TREATMENT.**

CW	Type	P/N	Box (Units)	Patient Size/ Weight	Wrap Height/Width (m)
Cure Wrap	Adult	508-03500	Box (X8)	168-180cm (over 66")	2.030/1.354
	Infant	508-03518	Box (X8)	Up to 4Kg	0.659/0.448
		508-03521	Box (X8)	4-7Kg	0.698/0.602
Pediatric Cure Wrap	Small	PED-SM008	Box (x8; multi size)	Up to 4Kg (X4), 4-7Kg (X4)	0.659/0.448 0.698/0.602
	Medium	PED-MD008	Box (x8; multi size)	7-11 Kg (X4), 79-91cm (X4)	0.981/0.628 1.118/0.740
	Large	PED-LA008	Box (x8; multi size)	91-104cm (X4), 104-122 cm (X4)	1.225/ 0.841 1.390/ 1.054
	X-Large	PED-XL008	Box (x8; multi size)	122-135cm (X4), Over 135cm (X4)	1.582/1.1193 2.030/1.354

## One-Piece

Body shaped garment  
Easy Application  
Facilitates wrapping of your patient

## 85% body area coverage

Optimal conductive energy transfer

## Maximum heat transfer

High velocity flow rate for maximum heat transfer

## Precision

ThermoRegulates in concert with CritiCool

## Soft mesh garment

Biocompatible, Latex Free

## Uncomplicated

Secured by Hook and Loop strips

## Re-apply

Permits patient handling and skin supervision.

## Safe

Irritation free

## Anywhere, Anytime

Safe in X-ray, Cath Lab, CT scan, MRI and during defibrillation

## Risk Free

Avoid risk of infections and complications associated with invasive cooling options

## Clutter free

One tube (with Double or Triple inlets)

## Disposable

No infection issues

## We fit all

- Varying Pediatric and Neonatal sizes, according to Weight and Height  
- One size fits all adults



CE 0473 FDA Cleared