

Mediflex® FlexArm™ Plus

Stronger, Smoother, Easier To Use!

Mediflex®, innovator of the Original FlexArm™ Universal Surgical Holding & Positioning Arm, introduces its latest generation technology.

The New FlexArm™ Plus design now features:

- Sleek, rounded body and knobs
- Quicker, smoother tightening
- Greater Rigidity - up to 2 times stronger than previous designs



- **No more need to RESET FlexArm™ when performing routine cleaning & sterilization.**

The FlexArm™ Plus is designed for speed and strength. 3 handle turns makes the arm rigid and ready for holding – 3 turns counterclockwise prepares the arm for cleaning and sterilization.

The FlexArm™ Plus is available in many configurations - see reverse side.

Mediflex®
SURGICAL PRODUCTS

Mediflex® FlexArm™ to FlexArm™ Plus Cross Reference - Ordering Information

The New Mediflex® Universal FlexArm™ Plus can be ordered in the same system configurations as original Universal FlexArm™ Systems. Please refer to our catalog SHPA-0913, pages 4-7 for a complete description of all available systems and arm configurations.



FlexArm™	FlexArm
Ref. No.	Plus™
	Ref. No.

Standard Single Arm Systems Description

69054-QC	99054-QC	Single Arm System with Quick Connect Hex Fitting
69054	99054	Single Arm System with Fixed Stainless Steel Tip
69054-D	99054-D	Single Arm System with Fixed Delrin Tip

Standard Twin Arm Systems

69055-QC	99055-QC	Twin Arm System with Quick Connect Hex Fittings
69055	99055	Twin Arm System with Fixed Stainless Steel Tips
69055-D	99055-D	Twin Arm System, with Fixed Delrin Tips
69055-DS	99055-DS	Twin Arm System, with One Fixed Stainless Steel Tip & One Fixed Delrin Tip

Long Reach Single Arm Systems

69054-QCLR	99054-QCLR	Single Arm System, Long Reach Post with Quick Connect Hex Fitting
69054-LR	99054-LR	Single Arm System, Long Reach with Fixed Stainless Steel Tip
69054-DLR	99054-DLR	Single Arm System, Long Reach with Fixed Delrin Tip

Long Reach Twin Arm Systems

69055-QCLR	99055-QCLR	Twin Arm, Long Reach Post with 2 Quick Connect, Hex Fittings
69055-LR	99055-LR	Twin Arm System, Long Reach, with Fixed Stainless Steel Tips
69055-DLR	99055-DLR	Twin Arm System, Long Reach, with Fixed Delrin Tips
69055-DSLR	99055-DSLR	Twin Arm System, Long Reach, with One Fixed Stainless Steel Tip & One Fixed Delrin Tip

Companion FlexArm™/FlexArm Plus™

69056-QC	99056-QC	Arm Only with Quick Connect Hex Fitting
69056	99056	Arm Only with Permanent Stainless Steel Tip
69056-D	99056-D	Arm Only with Permanent Delrin Tip
69056-QC-BSO1	99056-QC-BSO1	Arm Only with Distal Swivel Option & Quick Connect Hex Fitting

For Additional Information Contact Mediflex
or your Local Mediflex Distributor



Mediflex®
SURGICAL PRODUCTS

www.mediflex.com

Toll-free: 1-800-764-0474

Phone: 1-631-582-6424

Fax: 1-631-582-4520

© 2015 Flexbar Machine Corporation. All rights reserved.