

HOLTERY EKG

1. CardiUp3 3 kanałowy
2. CardiUp3/12 3 lub 12 kan.
3. CardioMera do 5 kan.

CardiUP!12

3/12-kanałowy Holter ECG



CardiUP!

- możliwość rejestracji do 7 dni 3 lub 12 odprowadzeń
- duży wyświetlacz LCD
- przycisk zdarzeń
- programowanie ręczne lub z PC
- mini USB
- oprogramowanie profesjonalne wolne od zabezpieczeń
można instalować na wielu komputerach bez dodatkowych kosztów

CardiUP Holter ECG




CardioMera

podstawowy Holter ECG

CardioMera

- 1-, 2-, 3- lub 5-kanalowy zapis
- łatwa wymiana kabli pacjenta w razie uszkodzenia
- mały rozmiar i mała waga
- przycisk zdarzeń



Specifications	CardiUP!12	CardiUP!3	CardioMera
dimensions	67*93*21 mm		66*59*17 mm
weight (excl. batteries)	92g		50g
memory	internal		external (a card reader unit is needed for retrieving data)
data transfer	USB 2.0 (mini USB)		USB optical cable
power supply	1 AA battery		1 AAA battery
LCD display	48*28 mm		—
number of keys	4 (event, up, down, select)		1 (event)
number of channels	3/12 independent channels (7 or 10 leads)	3 independent channels (7 leads)	1, 2, 3, 5 channels (7 leads)
operation	up to 168 hours of continuous ECG		up to 96 hours of continuous ECG
sampling, A/D resolution	300, 600 or 1200 Hz, 12 bit		600 or 1200 Hz, 12 bit
storage quality	300, 600 or 1200 Hz		150, 300 or 600 Hz
conformity certificates	CE		

Meditech Ltd.
 Budapest 1184 Mikszáth Kálmán u. 24.
 Tel. 36(1)280-82-32 • meditech@meditech.eu



Holter ECG REPORT

PATIENT DATA

name	0001_0001	Meditech	date of birth	2001.01.01.
ID				

HOLTER SUMMARY

The measurement started	2018.03.09. 11:57	Total beats	96311	
Total measurement time	24 h 00 min	Normal	10553	(11,0%)
Valuable length:	24 h 00 min	VES	127	(0,1%)
		SVES	223	(0,2%)
		PM	85408	(88,7%)
Channels / sampling rate	12 / 300 Hz	Other	0	(0,0%)
Total lead off / noise total	00 h 00 min	Artifact	0	(0,0%)
Pacemaker spike detection	enabled			


16:46:48

A+

HR=82

HR

HRV TACHOGRAM



HRV BAR CHART



CardioVisions oprogramowanie – bezpłatne

Specyfikacje	CardiUP!2	CardiUP!3	CardioMera
accuracy	99,90%		
beat-by-beat ECG annotation	✓		
event distribution	✓		
ST analysis	✓		
QT analysis	✓		
PQ analysis	✓		
PM analysis	✓		
time domain HRV analysis	✓		
frequency domain HRV analysis	✓		
arrhythmia analysis	✓		
graphical and tabular display of data	✓		
customizable screens	✓		
customizable & printable report	✓		
personal & network use	✓		



Meditech Ltd • H-1184 Budapest, Mikszáth Kálmán u. 24.
phone: +36 1 280 8232, +36 1 280 8233 • email: meditech@meditech.eu
www.meditech.eu