



FRESH AIR UVC STERILIZE RESPIRATOR

Committed To Photodynamic Life Research
Protect Health With Light Biological Effect



PROTECT THE SAFE LIFE UNDER THE NORMALIZATION OF EPIDEMIC SITUATION

LAMPLIC Innovative Products

20 years of UV development experience

Caring for life and health



UVC TWO-WAY DISINFECTION ENSURE BREATHING MORE SAFER

Built-in Uvc Led Germicidal Light With
Photocatalyst Filteractively Kills Bacteria And
Viruses In Inhaled And Exhaled Air.



ACTIVE FAN AIR SUPPLY BREATHING COMFORTABLY AND FREELY

High-speed Mute Fan With Up To 15000
Revolutions Ensures That The Breathing
Space Is In A Positive Pressure State And
Breathing Is Comfortable And Free






INGENIOUS QUADRUPLE PROTECTION DESIGN CUT OFF THE SPREAD OF THE VIRUS

Photocatalyst / Activated Carbon Absorbs /
Nano Silver Ion / N95 Meltblown Cloth

Quadruple effect filter , fully ensure the sterilization effect

COOL AND STYLISH DESIGN LATEST POPULAR EQUIPMENT

Cool And Fashionable Design, Combining The Latest
UVC Germicidal Light Technology With Modern
Popular Elements, Perfect Combination Of
Technology And Fashion



EXTREME LIFE NEW WAY OF BREATHING

Medical Grade Silicone Mask And Ear Strap
Ensure The Perfect Fit Of The Face And The
Comfort And Safety Of The Skin.





LONG WORKING TIME CLEAN AIR SUPPLY

Committed To Photodynamic Life Research
Protect Health With Light Biological Effect

BREATHE FREELY MORE SPARE PARTS

Abundant Spare Parts Are Supplied To Meet
The Needs Of Cleaning, Replacement And Daily
Consumption, Making Safety More Secure.





Model	M100
Size	87x 112 x 66mm
Weight	145 g
Working Time	4 Hours
Wearing Range	480 ~ 570 mm
Efficiency	> N95
Executive Standard	Q31/0117000290C002-2017
Material	Silicone / Plastic
Type	Positive Pressure
Input	5V 500mA
Max Noise	< 50 dB
Warranty Time	1 year

