


# Bronchoscope Stabilization System



Providing secure, hands-free holding of flexible bronchoscopes during navigational bronchoscopy –while evaluating and managing airway disease

- Eliminates an assistant for scope holding
- Reduces radiation exposure
- Maintains scope position during imaging and intervention
- Helps to reduce trauma to the anatomy
- Shortens procedure time

## Bronchoscope Stabilization System

table-mounted articulating arm offers static, hands-free holding of flexible bronchoscopes. One handle controls all three joints for ease of rapid positioning or adjustment.

### Flexible Scope Tip

allows for rapid, secure attachment to flexible bronchoscopes. Innovative design incorporates replaceable Silicone Liners for scope stability and cleanability. (*Silicone Liners sold separately*)



Silicone Liners provide a soft contact for scope protection and increase grip strength for various scope configurations.

## ORDERING INFORMATION

38000	<b>Bronchoscope Stabilization System</b> <i>Includes Fixed Flexible Scope Tip and Table Rail Clamp</i>
73000-SA	<b>StrongArm™ Holding &amp; Positioning System</b> <i>Includes Quick-Connect Hex Fitting and Table Rail Clamp</i>
99706	<b>Flexible Scope Tip with Hex Fitting</b> <i>For attachment to Mediflex StrongArm™ Systems with QC Fitting</i>
070305	<b>Silicone Liners, 12 per pack, Reusable/Sterilizable</b>

**Mediflex® Surgical Products**  
250 Gibbs Road  
Islandia, NY 11749 USA

631.582.6424  
info@mediflex.com  
www.mediflex.com



Shown: 73000-SA  
with 99706.  
Scope not included.