




### VIDEO COLONOSCOPE **EC-530DM/DL** Therapeutic Treatment



These lower GI tract endoscopes have two working channels (3.8 mm / 2.8 mm), intended to be useful for treatments such as EMR.



Viewing direction	0° (Forward)
Field of view	140°
Observation range	3–100 mm
Bending capability	Up 180° /Down 180° Right 160° /Left 160°
Distal end diameter	12.8 mm
Flexible portion diameter	12.8 mm
Working channel diameter	3.8/2.8 mm
Working length	1,330/1,690 mm
Total length	1,645 mm/2,005 mm



### VIDEO SIGMOIDOSCOPE **ES-530WE**



ES-530WE is a sigmoidoscope with an effective length of 790 mm. The working channel diameter is 3.8 mm, and it is equipped with a water jet function.



Viewing direction	0° (Forward)
Field of view	140°
Observation range	3–100 mm
Bending capability	Up 180° /Down 180° Right 160° /Left 160°
Distal end diameter	12.8 mm
Flexible portion diameter	12.8 mm
Working channel diameter	3.8 mm
Working length	790 mm
Total length	1,090 mm

