

Disposable robotic tips 300 µl conductive

235903-H 300 µl



Product description

Conductive tips, 300 µl, filter, tray, blister-5, 96 pieces in a tray


Part Number	235903-H
GTIN	N.A.
GDMN	Pipet tip
Tip volume	300 µl
Color Tip	black (conductive)
Color Tray	clear
Sterile	no
Re-usable	no
Patient contact	no
Contact with blood and blood plasma	yes

Classification


Europe CE labeling	CE-IVD	<i>in accordance with Directive 98/79/EC on in vitro diagnostic medical devices as accessory article which, whilst not being an in vitro diagnostic medical device, is intended specifically by Ritter GmbH to be used together with a device to enable that device to be used in accordance with its intended purpose</i>
Purity		Produced under cleanroom conditions (ISO 14644-1 class 8) Free of RNase, DNase, DNA and PCR inhibitors.

Materials	
Tip	PP – Polypropylene
Tray	PP – Polypropylene
Filter	PE – Polyethylene

Tip Dimensions	
Tip Length	59.9 mm
Concentricity	≤ 1.6
Filter depth	9.0 mm



Tray Dimensions	
Tray length	122.0 mm
Tray width	82.6 mm



Packaging	
Blister	96 tips per tray, 5 trays per blister
Box	12 blister per box (5,760 tips)
Pallet	16 boxes (92,160 tips)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website
www.picov.eu

Manufactured in Germany by Ritter GmbH -Disposable robotic tips , 300µL , conductive
Specifications are subject to change without notice **CE** **IVD**