





Technical data

Performance data	P Version	M Version
Applications	Compact stacking of ARGUS infusion and syringe pumps, power supply to the pumps	Compact stacking of ARGUS infusion and syringe pumps, power supply to the pumps, connection to a central nurse call system, provide ability to connect pumps to hospital IT system, such Patient Data Management System (PDMS).
Cascading	NA .	Up to 2 docking stations
Pump connection slots	A100P: 5 slots A60P: 3 slots	A100M: 5 slots A60M: 3 slots
Pump compatibility	A600S, A606S, A70XV, A71XV	A600S, A606S, A70XV, A71XV
History and Interface		
History and Event Log Nurse call	NA NA	200 events saved in real time Relay 24 V / 0.2 A, configurable pulse, potential-free
RS-232 interface	NA NA	Electrically isolated, connectable to PDMS
Safety		
Alarms Pump Monitoring	NA NA	Accoustic and visual Configurable that nurse call is activated when inserted pump has an alarm
Power supply		
Mains power supply Power consumption	230 VAC ±10% (option:115 VAC ±10%), 50–60 Hz 50 VA maximum	230 VAC ±10% (option:115 VAC ±10%), 50–60 Hz 50 VA maximum
Product Specification		
Dimensions	A100P: 84 x 602 x 80 mm (W x H x D) A60P: 84 x 420 x 80 mm (W x H x D)	A100M: 84 x 602 x 80 mm (W x H x D) A60M : 84 x 420 x 80 mm (W x H x D)
Housing Operating temperature	PC/ABS high performance plastic +5°C to +40 °C	PC/ABS high performance plastic +5°C to +40 °C
Storage temperature	±0°C to +40°C	±0°C to +40°C
Transport temperature	-10°C to +50°C 20–90 % without condensation	-10°C to +50°C 20–90 % without condensation
Relative humidity Weight	A100P: 2.3 kg	A100M: 2.4 kg
3	A60P: 1.9 kg	A60M: 2.0 kg
Maintenance	Safety inspection every 24 months	Safety inspection every 24 months
Compliance		
IP protection class Protection class	IPX1 Class I(level of protection against electric shock)	IPX1 Class I(level of protection against electric shock)
MDD classification		llb
Regulation and Electrical Safety	EN 60601-1	EN 60601-1, EN 60601-1-4
Electromagnetic compatibility Certification CODAN ARGUS AG	EN 60601-1-2 ISO 13485, ISO 9001	EN 60601-1-2 ISO 13485, ISO 9001
Certification device	CE 0120, directive 93/42/EEC	CE 0120, directive 93/42/EEC
Accessories		
Attachment system Software	2 combination clips, suitable for pole- and rail systems NA	2 combination clips, suitable for pole and rail systems ARGUSservice: program for simple maintenance and configu
SULWATE	IVA	ration
Cables	NA	Cable for PC / PDMS / cascading

Possible combinations of docking stations and ARGUS pumps





CODAN ARGUS AG, Oberneuhofstrasse 10, CH-6304 Baar Phone +41 41 785 09 44, Fax +41 41 785 09 40 codan@codanargus.com, www.codanargus.com