

# SWISS LITHOPUMP®

→ FOR ALL LITHOCLAST® DEVICES



**VENTURI  
PRINCIPLE**

# PRECISELY ADJUSTABLE → FOR ALL SUCTION REQUIREMENTS OF LITHOTRIPSY PROCEDURES

## FAILSAFE AND ULTRAQUIET

- > Runs on compressed air (venturi principle, no electricity)

## EASE OF USE

- > Plugs into the compressed air connection of the Swiss LithoClas<sup>®</sup>

## CLEAN AND HYGENIC

- > Double bacteria retention filter for increased hygiene safety
- > Minimum air consumption due to pneumatic auto control circuit

### SWISS LITHOPUMP<sup>®</sup>

Swiss LithoPump<sup>®</sup> basic unit including compressed air connection, filters and connection tubing

**FT-184#A**

- > Swiss LithoPump<sup>®</sup> with Swiss LithoClas<sup>®</sup> Master



### ACCESSORIES

Overflow protection filters (10x) **DT-096**

Disposable suction bags 2l (25x) **DT-097**

Double reusable suction jar set 2l **EP-035**

Autoclavable suction jar **EP-036**

Reusable suction jar **EP-037**

(for disposable suction bags)

Clamp holder **EP-038**

Change-over valve **EG-088**

Cart **DW-031**



EMS Electro Medical Systems S.A.  
Chemin de la Vuarpillièrè 31  
1260 Nyon  
Switzerland

Tel. +41 22 99 44 700  
Fax +41 22 99 44 701  
welcome@ems-ch.com  
www.ems-medical.com

**EMS**<sup>+</sup>  
ELECTRO MEDICAL SYSTEMS